


EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUSUBAULT AGENCY GOODYEAR PD
DATE 8/5/2010 TIME 1042
INTOXILYZER SERIAL # 80-000317 LOCATION GOODYEAR PD VENEZIA

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 101 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE _____


1144

INTOXILYZER 8000

Location: GOODYEAR PD VENIDA

Serial Number: 80-000312

Core Version: 8105.48

08/05/2010

10:42:25

Standard Lot#: 672733

Last Changed By: M. GUILBAULT #1144

CAS: M. GUILBAULT #1144

GOODYEAR PD

Operator: M. GUILBAULT #1144

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	10:43:36
0.100 Cal Check	0.101	10:43:57
Air Blank	0.000	10:44:27

Preventative Maintenance

Performed -- Timer Reset