

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUSCIBAUT 1144 AGENCY Goodyear PD
DATE 2/10/2010 TIME 2306
INTOXILYZER SERIAL # 80-000312 LOCATION GPD

- (C) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (S) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(S) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
(S) 4. Air blank completed.
(S) 5. Calibration check completed. Test results 0. 101 AC.
(S) 6. Air blank completed.
(S) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE _____
M. Guscibaut 1144

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-000312
02/10/2010

23:06:35

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	23:07:02
Cal Check	0.101	23:07:21
Air Blank	0.000	23:07:51


Operator's Signature

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GURUBAULT 1144 AGENCY GOODYEAR PD
DATE 2/24/2010 TIME 1739
INTOXILYZER SERIAL # 80-00313 LOCATION STATION #2

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 101 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE _____

M. Gurubault 1144


INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-000312
02/24/2010

17:39:20

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	17:39:47
Cal Check	0.101	17:40:06
Air Blank	0.000	17:40:36



Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. GUSUBAULT 1144 AGENCY GOODYEAR PD
DATE 2/24/2010 TIME 1740
INTOXILYZER SERIAL # 80-000 3/2 LOCATION STATION #2

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.


OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

AMBIENT FADL DURING AIR BLANK / SPRAYED
MOUTH WAS IN THE AREA.

SIGNATURE  1144

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number : 80-000312
02/24/2010

17:40:42

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST
DOB: 10/10/2000
Sex: M Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	17:41:31
Diagnostic Test	Pass	17:41:59
Air Blank	0.000	17:42:26
0.100 Cal Check	0.102	17:42:45
Air Blank	0.000	17:43:14
Subject Test	0.000	17:43:30
Air Blank	0.000	17:43:57
Five Minute Wait	ABT*	
Air Blank	0.000	17:44:31

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number : 80-000312
02/24/2010

17:44:35

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST
DOB: 10/10/2000
Sex: M Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	17:45:08
Air Blank	RFI**	17:45:21

*Improper Sample
**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number : 80-000312
02/24/2010

17:45:44

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST
DOB: 10/10/2000
Sex: M Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	17:46:34
Diagnostic Test	Pass	17:47:02
Air Blank	AMB*	17:47:30
Air Blank	0.000	17:48:01

*Ambient Fail

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number : 80-000312
02/24/2010

17:48:20

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST
DOB: 10/10/2000
Sex: M Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	17:49:08
Diagnostic Test	Pass	17:49:35
Air Blank	0.000	17:50:02
0.100 Cal Check	0.099	17:50:21
Air Blank	0.000	17:50:51
Subject Test	0.089*	17:53:59
Air Blank	0.000	17:54:31
Air Blank	0.000	17:54:58
Subject Test	XXX**	17:55:27
Air Blank	0.000	17:56:01

*Deficient Sample
**Invalid Sample
Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUSLBOUT 1144 AGENCY GOODYEAR PD
DATE 2/24/2010 TIME 1:1756
INTOXILYZER SERIAL # 80-000512 LOCATION STATION #2

- ① 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- ~~2.~~ 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
~~3.~~ 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
~~4.~~ 4. Air blank completed.
~~5.~~ 5. Calibration check completed. Test results 0. 099 AC.
~~6.~~ 6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

INTOXILYZER 8000

Location: GOODYEAR PD

Serial Number: 80-000312

02/24/2010

17:56:23

Standard Lot#: 672733

Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	17:56:51
Cal Check	0.099	17:57:10
Air Blank	0.000	17:57:39



Operator's Signature