

EXHIBIT 2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

AGENCY GOODYEAR POLICE DEPARTMENT
DATE _____ TIME 1803
INTOXILYZER SERIAL NO. 80-000312 LOCATION GOODYEAR PD
QA SPECIALIST Will
(print name)

- (14) 1. Ensure dry gas tank is attached to instrument and contains a known alcohol standard 100 AC.

OR

Pour a standard alcohol solution of known value, _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (2) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (2) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (2) 4. Air blank completed.
- (2) 5. Calibration check completed. Test results 0. 099 AC.
- (2) 6. Air blank completed.
- (2) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE: _____



INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number: 80-000312

01/21/2006

18:03:45

Standard Lot#: 431001G

Last Changed By: W. WILL #974

QAS: W. WILL #974

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	18:04:13
Cal Check	0.099	18:04:31
Air Blank	0.000	18:05:01

Operator's Signature

EXHIBIT 3

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

AGENCY GOODYEAR POLICE DEPARTMENT
DATE 1/21/06 TIME 1740
INTOXILYZER SERIAL NO. 80-000312 LOCATION GOODYEAR PD
QA SPECIALIST WILL
(Print Name)

(2) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

(2) 1. Clock time check.
(2) 2. Date check.

OPERATIONAL TESTS

(2) 1. Alcohol-free subject test result 0. 000 AC.
(2) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
(2) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
(2) 4. Calibration standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES X NO _____

COMMENTS: _____

SIGNATURE _____

EXHIBIT 2

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE


AGENCY GOODYEAR POLICE DEPARTMENT
DATE 1/21/2006 TIME 17:38
INTOXILYZER SERIAL NO. 80-000712 LOCATION GOODYEAR PD
QA SPECIALIST Will
(print name)

- (X) 1. Ensure dry gas tank is attached to instrument and contains a known alcohol standard 100 AC.

OR

Pour a standard alcohol solution of known value, _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 100 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE: _____



INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number: 80-000312

01/21/2006

17:38:58

Standard Lot#: 431001G

Last Changed By: W. WILL #974

QAS: W. WILL #974

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	17:39:26
Cal Check	0.100	17:39:44
Air Blank	0.000	17:40:14

Operator's Signature

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-000312

01/21/2006

17:40:39

Standard Lot#: 431001G

Last Changed By: W. WILL

#974

QAS: W. WILL

#974

GOODYEAR PD

Operator: W. WILL

#974

GOODYEAR PD

Subject: TEST,TEST

DOB: 11/11/1111

Sex: M

Weight: 1

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	17:41:32
Diagnostic Test	Pass	17:42:00
Air Blank	0.000	17:42:27
0.100 Cal Check	0.099	17:42:46
Air Blank	0.000	17:43:15
Subject Test	0.000	17:43:29
Air Blank	0.000	17:43:57
Five Minute Wait	ABT*	
Air Blank	0.000	17:44:44

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-000312

01/21/2006

17:45:42

Standard Lot#: 431001G

Last Changed By: W. WILL

#974

QAS: W. WILL

#974

GOODYEAR PD

Operator: W. WILL

#974

GOODYEAR PD

Subject: TEST,TEST

DOB: 11/11/1111

Sex: M

Weight: 1

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	17:46:17
Air Blank	0.000	17:46:46

*Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-000312

01/21/2006

17:46:52

Standard Lot#: 431001G

Last Changed By: W. WILL

#974

QAS: W. WILL

#974

GOODYEAR PD

Operator: W. WILL

#974

GOODYEAR PD

Subject: TEST,TEST

DOB: 11/11/1111

Sex: M

Weight: 1

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	17:47:35
Air Blank	0.000	17:48:04

*RFI Detect

Not a Successfully Completed Test Sequence

EXHIBIT 3

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R9-14-404(A)

ARIZONA DEPARTMENT OF HEALTH SERVICES

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

AGENCY GOODYEAR POLICE DEPARTMENT

DATE _____ TIME _____

INTOXILYZER SERIAL NO. 80-000312 LOCATION GOODYEAR PD

QA SPECIALIST WILL

(Print Name)

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

() 1. Clock time check.

() 2. Date check. _____

OPERATIONAL TESTS

() 1. Alcohol-free subject test result 0. _____ AC.

() 2. Error recognition logic system functioning.

Not a Successfully Completed Test Sequence printed _____

() 3. Proper sample recognition system.

Not a Successfully Completed Test Sequence printed _____

Deficient sample printed.

() 4. Calibration standard 0. _____ AC. Result 0. _____ AC.

Instrument is operating properly and accurately. YES _____ NO _____

COMMENTS: _____

SIGNATURE _____

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-000312

01/21/2006

17:55:09

Standard Lot#: 431001G

Last Changed By: W. WILL

#974

QAS: W. WILL

#974

GOODYEAR PD

Operator: W. WILL

#974

GOODYEAR PD

Subject: TEST,TEST

DOB: 11/11/1111

Sex: M

Weight: 1

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	17:56:18
Diagnostic Test	Pass	17:56:45
Air Blank	0.000	17:57:12
0.100 Cal Check	0.099	17:57:31
Air Blank	0.000	17:58:01
Subject Test	0.090*	18:01:09
Air Blank	0.000	18:01:42
Air Blank	0.000	18:02:10
Subject Test	XXX**	18:02:38
Air Blank	0.000	18:03:12

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-000312

01/21/2006

17:48:09

Standard Lot#: 431001G

Last Changed By: W. WILL

#974

QAS: W. WILL

#974

GOODYEAR PD

Operator: W. WILL

#974

GOODYEAR PD

Subject: TEST,TEST

DOB: 11/11/1111

Sex: M

Weight: 1

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	17:48:57
Diagnostic Test	Pass	17:49:24
Air Blank	0.000	17:49:52
0.100 Cal Check	0.100	17:50:10
Air Blank	0.000	17:50:40
Subject Test	0.049*	17:53:48
Air Blank	RFI**	17:54:16
Air Blank	0.000	17:54:45

*Deficient Sample

**RFI Detect

*TURNED ON CELLPHONE
WHEN UNIT WAS AIR BLANKING*

Not a Successfully Completed Test Sequence