

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. ROGERS AGENCY GOODYEAR P.D.
DATE 6.20.07 TIME 1854
INTOXILYZER SERIAL # 80.000312 LOCATION GOODYEAR P.D.

- ~~FF~~ 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- ~~FF~~ 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- ~~FF~~ 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- ~~FF~~ 4. Air blank completed.
- ~~FF~~ 5. Calibration check completed. Test results 0. 099 AC.
- ~~FF~~ 6. Air blank completed.
- ~~FF~~ 7. Remove printed record. Attach the record to the completed checklist.

LOT 618601 I
EXP 7-7-08

SIGNATURE [Signature] #1123

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST J. ROGERS AGENCY GOODYEAR P.D
DATE 6-20-07 TIME 1856
INTOXILYZER SERIAL # 80-000312 LOCATION GOODYEAR PD

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

LOT 618601Z
EXP 7-7-08

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE [Signature] #1123

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-000312

06/20/2007

18:56:55

Standard Lot#: 618601I

Last Changed By: J. TRAINING

#99999

QAS: J. ROGERS

#1123

GOODYEAR PD

Operator: J. ROGERS

#1123

GOODYEAR PD

Subject: TEST,TEST

DOB: 01/01/0101

Sex: M

Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:57:45
Diagnostic Test	Pass	18:58:13
Air Blank	0.000	18:58:40
0.100 Cal Check	0.099	18:58:59
Air Blank	0.000	18:59:28
Subject Test	0.000	18:59:46
Air Blank	0.000	19:00:13
Five Minute Wait	ABT*	
Air Blank	0.000	19:00:50

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-000312

06/20/2007

19:00:58

Standard Lot#: 618601I

Last Changed By: J. TRAINING

#99999

QAS: J. ROGERS

#1123

GOODYEAR PD

Operator: J. ROGERS

#1123

GOODYEAR PD

Subject: TEST,TEST

DOB: 01/01/0101

Sex: M

Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	19:01:36
Air Blank	0.000	19:02:04

*Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-000312

06/20/2007

19:02:28

Standard Lot#: 618601I

Last Changed By: J. TRAINING

#99999

QAS: J. ROGERS

#1123

GOODYEAR PD

Operator: J. ROGERS

#1123

GOODYEAR PD

Subject: TEST,TEST

DOB: 01/01/0101

Sex: M

Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	19:03:07
Air Blank	0.000	19:03:37

*RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-000312

06/20/2007

19:04:06

Standard Lot#: 618601I

Last Changed By: J. TRAINING

#99999

QAS: J. ROGERS

#1123

GOODYEAR PD

Operator: J. ROGERS

#1123

GOODYEAR PD

Subject: TEST,TEST

DOB: 01/01/0101

Sex: M

Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	19:04:54
Diagnostic Test	Pass	19:05:22
Air Blank	0.000	19:05:49
0.100 Cal Check	0.099	19:06:08
Air Blank	0.000	19:06:37
Subject Test	0.016*	19:09:46
Air Blank	0.000	19:10:16
Air Blank	0.000	19:10:43
Subject Test	XXX**	19:11:10
Air Blank	0.000	19:11:43

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. ROGERS AGENCY GOODYEAR P.D.
DATE 6-20-07 TIME 1912
INTOXILYZER SERIAL # 80-000312 LOCATION GOODYEAR P.D.

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0. 099 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

LOT 618601 I
EXP 7-7-08

SIGNATURE [Signature] #1123

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number: 80-000312

06/20/2007

19:12:29


Standard Lot#: 618601I

Last Changed By: J. TRAINING #99999

QAS: J. ROGERS #1123

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	19:12:57
Cal Check	0.099	19:13:15
Air Blank	0.000	19:13:45



Operator's Signature