

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GLEIBAUER 1144 AGENCY GREENE PA
DATE 3/13/2010 TIME 1533
INTOXILYZER SERIAL # 80-000312 LOCATION GAS #2

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution. 100 AC
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$.
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test result 0.101 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-000312
03/13/2018

15:33:23

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	15:33:50
Cal Check	0.101	15:34:09
Air Blank	0.000	15:34:38


Operator's Signature

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GIBSON AGENCY POUNCEY PD
DATE 3/18/2010 TIME 1838
INTOXILYZER SERIAL # 80000312 LOCATION GAD #2

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0.101 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-000312
03/18/2010

18:38:34

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	18:39:01
Cal Check	0.101	18:39:20
Air Blank	0.000	18:39:49



Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. G. ZAVIT 1144 AGENCY Goodyear PD
DATE 3/18/2010 TIME 1840
INTOXILYZER SERIAL # 80-000312 LOCATION CPD #2

Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 101 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

MOUTH ALCOHOL - DID NOT WAIT LONG ENOUGH BEFORE GIVING
SAMPLE.

SIGNATURE M. G. ZAVIT 1144

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number : 80-000312
03/18/2010 18:40:22

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST
DOB: 01/01/2000
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:41:10
Diagnostic Test	Pass	18:41:38
Air Blank	0.000	18:42:05
0.100 Cal Check	0.101	18:42:24
Air Blank	0.000	18:42:53
Subject Test	0.000	18:43:17
Air Blank	0.000	18:43:44
Five Minute Wait	ABT*	
Air Blank	0.000	18:44:19

*Sequence Aborted



Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: GOODYEAR PD

Serial Number : 80-000312

03/18/2010

18:44:45

Standard Lot#: 672733

Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144

GOODYEAR PD

Operator: M. GUILBAULT #1144

GOODYEAR PD

Subject: TEST,TEST

DOB: 01/01/2000

Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	18:45:22
Air Blank	RFI**	18:45:31

*Improper Sample

**RFI Detect

HA

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number : 80-000312
03/18/2010

18:45:57

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST
DOB: 01/01/2000
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:46:44
Diagnostic Test	Pass	18:47:12
Air Blank	0.000	18:47:39
0.100 Cal Check	0.101	18:47:58
Air Blank	0.000	18:48:28
Subject Test	0.111*	18:51:36
Air Blank	0.000	18:52:08
Air Blank	0.000	18:52:35
Subject Test	0.272	18:53:08
Air Blank	0.000	18:53:43
Five Minute Wait	ABT**	
Air Blank	0.000	18:54:43

*Deficient Sample
**Sequence Aborted

H

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number : 80-000312
03/18/2010 18:55:07

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST, TEST
DOB: 01/01/2000
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:55:52
Diagnostic Test	Pass	18:56:20
Air Blank	0.000	18:56:47
0.100 Cal Check	0.101	18:57:06
Air Blank	0.000	18:57:36
Subject Test	0.051	18:58:17
Air Blank	ABT*	18:58:31
Air Blank	0.000	18:58:59

*Sequence Aborted

fb

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number : 80-000312
03/18/2010

18:59:36

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST
DOB: 01/01/2000
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	19:00:25
Diagnostic Test	Pass	19:00:53
Air Blank	0.000	19:01:20
0.100 Cal Check	0.101	19:01:39
Air Blank	0.000	19:02:09
Subject Test	0.063	19:02:38
Air Blank	ABT*	19:02:47
Air Blank	0.000	19:03:15

*Sequence Aborted



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number : 80-000312
03/18/2010

19:03:21

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST, TEST
DOB: 01/01/2000
Sex: M

Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	19:04:12
Diagnostic Test	Pass	19:04:39
Air Blank	0.000	19:05:07
0.100 Cal Check	0.100	19:05:25
Air Blank	0.000	19:05:55
Subject Test	0.538	19:06:13
Air Blank	0.000	19:06:52
Five Minute Wait	ABT*	
Air Blank	0.000	19:07:24

*Sequence Aborted



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number : 80-000312
03/18/2010 19:07:29

Standard Lot#: 672733
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST, TEST
DOB: 01/01/2000
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	19:08:15
Diagnostic Test	Pass	19:08:42
Air Blank	0.000	19:09:09
0.100 Cal Check	0.101	19:09:28
Air Blank	0.000	19:09:58
Subject Test	XXX*	19:10:15
Air Blank	0.000	19:10:50

*Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GIEBELAULT AGENCY GOODYFAR PD
DATE 3/18/2010 TIME 1911
INTOXILYZER SERIAL # Y0-000312 LOCATION CPD #2

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 100 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

DPS Form Exh G-2 (Rev 05-01)

INTOXILYZER 8000

Location: GOODYEAR PD

Serial Number: 80-000312

03/18/2010

19:11:33

Standard Lot#: 672733

Last Changed By: M. GUILBAULT #1144

OAS: M. GUILBAULT #1144

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	19:12:00
Cal Check	0.100	19:12:19
Air Blank	0.000	19:12:49



Operator's Signature

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST N.W. AGENCY GOOD/EAR PD
DATE 3-19-2010 TIME 1036
INTOXILYZER SERIAL # 80000312 LOCATION GOOD/EAR PD

- (W) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (W) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (W) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (W) 4. Air blank completed.
- (W) 5. Calibration check completed. Test results 0. 1.102 AC.
- (W) 6. Air blank completed.
- (W) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000
Instrument Initialization
10:28:31 03/19/2010

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-000312
Core Version: 8105.48
03/19/2010 10:36:53

Standard Lot#: 672733
Last Changed By: W. WILL #974

CAS: W. WILL #974
GOODYEAR PD

Operator: W. WILL #974
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	10:37:26
0.100 Cal Check	0.102	10:37:47
Air Blank	0.000	10:38:17

Preventative Maintenance
Performed -- Timer Reset

Lab Use Only	
Downloaded?	Yes No
Date Updated	
Spreadsheet?	Inst. QAS

Intoxilyzer 8000 Update Form

Instrument Serial # 80-000312

Permanent Location GOODYEAR PD

COBRA Downloaded Not Downloaded

Primary QAS Information

Name N. WILL Badge 974

Agency GOODYEAR PD

Agency Address PO Box 5100
GOODYEAR AZ 85338

Work Phone 623 882 7093 Cell Phone 623 369 30072

Email nwill@GOODYEARAZ.GOV

Secondary QAS Information

Name M. GUILBAULT Badge 1144

Agency _____

Agency Address _____

Work Phone _____ Cell Phone _____

Email M.GUILBAULT@GOODYEARAZ.GOV

* If this instrument is taken care of by additional QAS officers, please put their information on the back of this form.

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. Lewis 1144 AGENCY GOODYEN PA
DATE 3/27/2010 TIME 0307
INTOXILYZER SERIAL # 80-000312 LOCATION 6A0 #2

- (H) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (S) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (S) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (H) 4. Air blank completed.
- (S) 5. Calibration check completed. Test results 0. 101 AC.
- (S) 6. Air blank completed.
- (S) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

INTOXLYZER 8300
Location: GOODYEAR PD
Serial Number: 80-000312
Core Version: 8135.48
03/27/2016

03:07:55

Standard Lot#: 672733
Last Changed By: D. WILK #974

DAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	03:08:25
0.100 Cal Check	0.100	03:09:46
Air Blank	0.000	03:05:15

Preventative Maintenance
Performed -- Timer Reset

