

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUILLMANT 1144 AGENCY GOODYEAR PD  
DATE 3/14/2011 TIME 0735  
INTOXILYZER SERIAL # 80-000312 LOCATION 10701 S ESTRELLA PRVY

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 102 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE \_\_\_\_\_ 1144

INTOXILYZER 8000  
Location: 10701 S ESTRELLA  
Serial Number: 80-000312  
Core Version: 8105.48  
03/14/2011

07:35:29

-----  
Standard Lot#: 16610100A1  
Last Changed By: M. GUILBAULT #1144  
-----

QAS: M. GUILBAULT #1144  
GOODYEAR PD  
-----

Operator: M. GUILBAULT #1144  
GOODYEAR PD  
-----

Test	g/210L	Time
Air Blank	0.000	07:36:09
0.100 Cal Check	0.102	07:36:30
Air Blank	0.000	07:36:59

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. GURBAUT 1144 AGENCY GOODYEAR PD  
DATE 3/14/2011 TIME 0737  
INTOXILYZER SERIAL # 80-000312 LOCATION 10701 SESTREUM PRING

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.  
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.  
 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 101 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

---

---

---

---

SIGNATURE M. Gurbaut 1144

INTOXILYZER 8000  
Location: 10701 S ESTRELLA  
Serial Number: 80-000312  
Core Version: 8105.48  
03/14/2011

07:37:21

Standard Lot#: 16610100A1  
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Subject: TEST,TEST  
DOB: 01/01/1990  
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	07:38:11
Diagnostic Test	Pass	07:38:46
Air Blank	0.000	07:39:13
0.100 Cal Check	0.103	07:39:34
Air Blank	0.000	07:40:04
Subject Test	0.000	07:40:16
Air Blank	0.000	07:40:43
Wait	ABT*	
Air Blank	0.000	07:41:17

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: 10701 S ESTRELLA  
Serial Number: 80-000312  
Core Version: 8105.48  
03/14/2011

07:41:21

---

Standard Lot#: 16610100A1  
Last Changed By: M. GUILBAULT #1144

---

QAS: M. GUILBAULT #1144  
GOODYEAR PD

---

Operator: M. GUILBAULT #1144  
GOODYEAR PD

---

Subject: TEST,TEST  
DOB: 01/01/1990  
Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	07:41:54
Air Blank	RFI**	07:42:02

\*Improper Sample  
\*\*RFI Detect

*[Handwritten signature]*

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: 10701 S ESTRELLA  
Serial Number: 80-000312  
Core Version: 8105.48  
03/14/2011

07:42:11

Standard Lot#: 16610100A1  
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Subject: TEST,TEST  
DOB: 01/01/1990  
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	07:42:58
Diagnostic Test	Pass	07:43:33
Air Blank	0.000	07:44:01
0.100 Cal Check	0.102	07:44:22
Air Blank	0.000	07:44:51
Subject Test	0.166*	07:47:59
Air Blank	0.000	07:48:31
Air Blank	0.000	07:48:58
Subject Test	XXX**	07:49:26
Air Blank	0.000	07:50:00

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUBELMANT 1144 AGENCY GOODYEAR PD  
DATE 3/14/2011 TIME 0750  
INTOXILYZER SERIAL # 80-000312 LOCATION 10701 SE STELLA PKWY

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 101 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE \_\_\_\_\_

INTOXILYZER 8000  
Location: 10701 S ESTRELLA  
Serial Number: 80-000312  
Core Version: 8105.48  
03/14/2011

07:50:14

-----  
Standard Lot#: 16610100A1  
Last Changed By: M. GUILBAULT #1144  
-----

QAS: M. GUILBAULT #1144  
GOODYEAR PD  
-----

Operator: M. GUILBAULT #1144  
GOODYEAR PD  
-----

Test	g/210L	Time
Air Blank	0.000	07:50:43
0.100 Cal Check	0.101	07:51:04
Air Blank	0.000	07:51:34

Preventative Maintenance  
Performed -- Timer Reset