

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-18-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUILBAULT AGENCY GOODYEAR PD
DATE 5/26/2010 TIME 1620
INTOXILYZER SERIAL # 80-000318 LOCATION GOODYEAR POLICE

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed, Test results 0. 101 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

ISTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-000312
Core Version: 8105.13
05/26/2010

16:25:07

Standard Lot#: 672733
Last Changed By: W. WILL 1974

QAS: M. CULBRAULT 41144
GOODYEAR PD

Operator: M. CULBRAULT 41144
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	16:25:46
0.100 Cal Check	0.101	16:26:10
Air Blank	0.000	16:26:20

Preventative Maintenance
Performed -- Timer Reset



INTOX 80-000312

**CERTIFICATE OF ANALYSIS
EBS - ETHANOL BREATH STANDARD**

CMI
316 EAST NINTH STREET
OWENSBORO, KY 42303

INVOICE#: 37420686
PO#: 374475
CMI ITEM #: 340135
DATE: 26 May 2010

ANALYTICAL ACCURACY: +/-0.002 BrAC or +/-2% whichever is greater.
PRODUCT EXPIRATION: 21 May 2012

CALGAZ LOT#: 767153
ETHANOL IN NITROGEN

COMPONENT	PPM	(BrAC)
ETHANOL	260.5	(0.100)
NITROGEN	BAL	

REFERENCE STANDARD	CYLINDER	CONCENTRATION PPM
NIST TRACEABLE STANDARDS*	D247819	259.6

* CERTIFICATION TRACEABLE TO N.I.S.T. RGM ETHANOL STANDARDS

MANUFACTURED DATE: 21 May 2010

CALGAZ CYLINDER SIZE: 6D

APPROVED BY : 

"We certify that all the cylinders for the Lot numbers identified herein are manufactured and tested within the requirements of CFR 49 part 178.65 and that physical and chemical test reports are on file and copies will be furnished upon request."

CALGAZ, Div. of Air Liquide Advanced Technologies U.S. LLC
821 Chesapeake Drive, Cambridge, MD 21613-0149
Phone: (410)228-6400 Fax: (410)228-4251

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. Gervault 1144 AGENCY GOODYEAR PD
DATE 5/31/2010 TIME 1525
INTOXILYZER SERIAL # 80-000312 LOCATION LAKE PLEASANT COMMAND VAN

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0. 100 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

INTOXILYZER 8000

Location: GOODYEAR POLICE VAN

Serial Number: 80-000312

Core Version: 8105-48

05/31/2010

15:23:23

Standard Lot#: 872733

Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144

GOODYEAR PD

Operator: M. GUILBAULT #1144

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	15:25:00
0.100 Cal Check	0.100	15:26:21
Air Blank	0.000	15:26:50

Preventative Maintenance
Performed -- Timer Reset