

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUZMAN AGENCY GOODYEAR PD
DATE 11/10/2009 TIME 0957
INTOXILYZER SERIAL # 80-000312 LOCATION GPD

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- ~~2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"~~
~~3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.~~
~~4. Air blank completed.~~
~~5. Calibration check completed. Test results 0. 105 AC.~~
~~6. Air blank completed.~~
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

INTOXILYZER 8000

Location: GOODYEAR PD

Serial Number: 80-000312

11/10/2009

09:56:44

Standard Lot#: 562576

Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	09:57:11
Cal Check	0.105 <i>th</i>	09:57:30
Air Blank	0.000	09:57:59



Operator's Signature