

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST W.T. RYAN AGENCY GOODYEAR P.D.
DATE 10-3-13 TIME 1559
INTOXILYZER SERIAL # 80-000312 LOCATION FIRE 2

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.10 AC.
Pour a standard alcohol concentration solution ^{OR} 0.20 AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0.098 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE W.T. Ryan 1100

INTOXILYZER 8000

Location: 10701 S ESTRELLA PKW

Serial Number: 80-000312

Core Version: 8105.48

10/03/2013

15:59:46

Standard Lot#: 32012100A1

Last Changed By: M. GUILBAULT #1144

QAS: W. RYAN #1100

GOODYEAR PD

Operator: W. RYAN #1100

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	16:00:18
0.100 Cal Check	0.098	16:00:38
Air Blank	0.000	16:01:08

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST W.T. RYAN AGENCY GOODYEAR P.D.
DATE 10-3-13 TIME 1603
INTOXILYZER SERIAL # 80-000312 LOCATION FIRE II

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 10 AC. Result 0. 98 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE W.T. Ryan 1100

INTOXILYZER 8000
Location: 10701 S ESTRELLA PKW
Serial Number: 80-000312
Core Version: 8105.48
10/03/2013 16:11:36

Standard Lot#: 32012100A1
Last Changed By: M. GUILBAULT #1144

QAS: W. RYAN #1100
GOODYEAR PD

Operator: W. RYAN #1100
GOODYEAR PD

Subject: TEST,TEST
DOB: 10/03/1959
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	16:12:58
Diagnostic Test	Pass	16:13:33
Air Blank	0.000	16:14:01
0.100 Cal Check	0.098	16:14:22
Air Blank	0.000	16:14:52
Subject Test	0.013*	16:17:59
Air Blank	0.000	16:18:29
Air Blank	0.000	16:18:56
Subject Test	XXX**	16:19:30
Air Blank	0.000	16:20:06

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: 10701 S ESTRELLA PKW
Serial Number: 80-000312
Core Version: 8105.48
10/03/2013 16:10:09

Standard Lot#: 32012100A1
Last Changed By: M. GUILBAULT #1144

QAS: W. RYAN #1100
GOODYEAR PD

Operator: W. RYAN #1100
GOODYEAR PD

Subject: TEST,TEST
DOB: 10/03/1959
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	16:10:51
Air Blank	0.000	16:11:19

*RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: 10701 S ESTRELLA PKW
Serial Number: 80-000312
Core Version: 8105.48
10/03/2013 16:03:02

Standard Lot#: 32012100A1
Last Changed By: M. GUILBAULT #1144

QAS: W. RYAN #1100
GOODYEAR PD

Operator: W. RYAN #1100
GOODYEAR PD

Subject: TEST,TEST
DOB: 10/03/1959
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	16:04:20
Diagnostic Test	Pass	16:04:54
Air Blank	0.000	16:05:22
0.100 Cal Check	0.099	16:05:43
Air Blank	0.000	16:06:12
Subject Test	0.000	16:06:29
Air Blank	0.000	16:06:57
Wait	ABT*	
Air Blank	0.000	16:07:46

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: 10701 S ESTRELLA PKW
Serial Number: 80-000312
Core Version: 8105.48
10/03/2013 16:08:24

Standard Lot#: 32012100A1
Last Changed By: M. GUILBAULT #1144

QAS: W. RYAN #1100
GOODYEAR PD

Operator: W. RYAN #1100
GOODYEAR PD

Subject: TEST,TEST
DOB: 10/03/1959
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	16:09:19
Air Blank	0.000	16:09:47

*Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: 10701 S ESTRELLA PKW

Serial Number: 80-000312

Core Version: 8105.48

10/03/2013

16:21:05

Standard Lot#: 32012100A1

Last Changed By: M. GUILBAULT #1144

QAS: W. RYAN #1100

GOODYEAR PD

Operator: W. RYAN #1100

GOODYEAR PD

Test g/210L Time

Air Blank 0.000 16:21:37
0.100 Cal Check 0.098 16:21:58
Air Blank 0.000 16:22:28

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST W.T. RYAN AGENCY GOODYEAR P.D.
DATE 10-16-13 TIME 1247
INTOXILYZER SERIAL # 80-000312 LOCATION STATION 4

(✓) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (✓) 1. Clock time check.
(✓) 2. Date check.

OPERATIONAL TESTS

- (✓) 1. Alcohol-free subject test result 0. 000 AC.
(✓) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
(✓) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
(✓) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE

W.T. Ryan 1100

INTOXILYZER 8000

Location: 10701 S ESTRELLA PKW

Serial Number: 66-000312

Core Version: 8165.46

10/15/2013

12:47:34

Standard Lot#: 32012100A1

Last Changed By: M. GUILSAULT #1144

DAS: W. RYAN #1100

GOODYEAR PD

Operator: W. RYAN #1100

GOODYEAR PD

Subject: TEST,TEST

DOB: 10/16/1980

Sex: M Weight: 180

15 Minute Depreciation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:48:27
Diagnostic Test	Pass	12:49:04
Air Blank	0.000	12:49:32
0.100 Cal Check	0.100	12:49:52
Air Blank	0.000	12:50:22
Subject Test	0.000	12:50:42
Air Blank	0.000	12:51:10
Wait	087*	
Air Blank	0.000	12:51:49

*Sequence Aborted

Not a Successfully
Completed
Test Sequence

INTOXILYZER 8000
Location: 19701 S ESTRELLA PKW
Serial Number: 80-003312
Core Version: 8195.48
10/16/2013 12:52:05

Standard Lot#: 32012100A1
Last Changed By: M. GUILBAULT #1144

QAS: W. RYAN #1100
GOODYEAR PD

Operator: W. RYAN #1100
GOODYEAR PD

Subject: TEST,TEST
DOB: 10/16/1980
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	12:52:55
Air Blank	0.000	12:53:23

*Improper Sample

Not a Successfully
Completed
Test Sequence

Not a Successfully Completed Test Sequence

MICROLYZER 8000
Location: 10701 S ESTABLER PKW
Serial Number: 80-000312
Core Version: 9105.48
10/16/2013
Standard Lot#: 3201210001
Last Changed By: M. GUILBERT #1144
Operator: M. RYAN #1100
Operator: M. RYAN #1106
GOODYEAR PD
DOB: 10/16/1980
Sex: M Weight: 180
Subject: TEST, TEST
15 Minute Depuration Period? Yes
Test Time: 9/21/08
Air Blank RFI* 0.000
Air Blank RFI* 0.000
12:54:37
12:55:35
*RFI Detect

INTOXILYZER 8000
Location: 10781 S ESTRELLA PKW
Serial Number: 00-000312
Core Version: 8105.48
10/16/2013 12:55:25

Standard Lot#: 32012100A1
Last Changed By: M. GUILBAULT #1144

DAS: U. RYAN #1100
GOODYEAR PD

Operator: U. RYAN #1100
GOODYEAR PD

Subject: TEST,TEST
DOB: 10/16/1980
Sex: M Weight: 180

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:56:35
Diagnostic Test	Pass	12:57:13
Air Blank	0.000	12:57:41
0.100 Cal Check	0.099	12:58:02
Air Blank	0.000	12:58:32
Subject Test	0.031*	13:01:39
Air Blank	0.000	13:02:10
Air Blank	0.000	13:02:35
Subject Test	XXX**	13:02:53
Air Blank	0.000	13:03:29

*Deficient Sample
**Invalid Sample
Begin new deprivation period

Not a Successfully
Completed
Test Sequence

INTOXILYZER 8000
Location: 10701 S ESTRELLA PKW
Serial Number: 80-000312
Core Version: 8105.48
10/16/2013

13:04:07

Standard Lot#: 32012100A1
Last Changed By: M. GUILBAULT #1144

DAS: W. RYAN #1100
GOODYEAR PD

Operator: W. RYAN #1100
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	13:04:39
0.100 Cal Check	0.098	13:05:00
Air Blank	0.000	13:05:30

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST W.T. RYAN 1100 AGENCY GOODYEAR P.D.
DATE 10-16-13 TIME 1244
INTOXILYZER SERIAL # 80-000312 LOCATION STATION 4

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 10 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 100 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE W.T. Ryan 1100

INDUSTRY 8000

Location: 10701 S ESTRELLA PKW

Serial Number: 80-800312

Core Version: 9105.48

10/16/2013

12:44:46

Standard Lot#: 32012:00A1

Last Changed By: M. GUILBALLT #1144

CAS: W. RYAN #1100

GOODYEAR PD

Operator: W. RYAN #1100

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	12:45:17
0.100 Cal Check	0.100	12:45:36
Air Blank	0.000	12:46:07

Preventative Maintenance
Performed -- Timer Reset