

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. ROGERS AGENCY GOODYEAR PD
DATE 9-1-06 TIME 1802
INTOXILYZER SERIAL # 80-000312 LOCATION GOODYEAR P.D.

~~1~~ 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution _____ AC.
OR

Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$

- ~~2~~ 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
~~3~~ 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
~~4~~ 4. Air blank completed.
~~5~~ 5. Calibration check completed. Test results 0. 099 AC.
~~6~~ 6. Air blank completed.
~~7~~ 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] #1123

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number: 80-000312

09/01/2006

18:02:58

Standard Lot#: 432401I

Last Changed By: J. VALLANDINGHAM #981

QAS: J. ROGERS #1123

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	18:03:26
Cal Check	0.099	18:03:44
Air Blank	0.000	18:04:14

 #1123

Operator's Signature

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. ROGERS AGENCY GOODYEAR P.D.
DATE 9-28-06 TIME 0555
INTOXILYZER SERIAL # 80-000312 LOCATION GOODYEAR P.D.

- ~~TS~~ 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- ~~TS~~ 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- ~~TS~~ 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- ~~TS~~ 4. Air blank completed. LOT 4324011
- ~~TS~~ 5. Calibration check completed. Test results 0. 099 AC. EXP 11/23/06
- ~~TS~~ 6. Air blank completed.
- ~~TS~~ 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] # 1123

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number: 80-000312

09/28/2006

05:55:10

Standard Lot#: 432401I

Last Changed By: J. VALLANDINGHAM #981

QAS: J. ROGERS #1123

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	05:55:38
Cal Check	0.099	05:55:56
Air Blank	0.000	05:56:26



Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST J. Rogers AGENCY GOODYEAR P.D.
DATE 9-28-06 TIME 0557
INTOXILYZER SERIAL # 80-000312 LOCATION GOODYEAR P.D.

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

LOT 4324011
EXP 11/23/06

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.

2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed


3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.

4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE



INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-000312

09/28/2006

05:57:30

Standard Lot#: 432401I

Last Changed By: J. VALLANDINGHAM #981

QAS: J. ROGERS #1123
GOODYEAR PD

Operator: J. ROGERS #1123
GOODYEAR PD

Subject: TEST,TEST


DOB: 01/01/0101

Sex: M Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	05:58:48
Diagnostic Test	Pass	05:59:16
Air Blank	0.000	05:59:43
0.100 Cal Check	0.099	06:00:02
Air Blank	0.000	06:00:31
Subject Test	0.000	06:00:51
Air Blank	0.000	06:01:18
Five Minute Wait	ABT*	
Air Blank	0.000	06:01:53

*Sequence Aborted

A handwritten signature in black ink, consisting of a stylized 'J' and 'R' intertwined, followed by the number '#1123' written to the right of the signature.

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-000312

09/28/2006

06:13:04

Standard Lot#: 432401I

Last Changed By: J. VALLANDINGHAM #981

QAS: J. ROGERS #1123
GOODYEAR PD

Operator: J. ROGERS #1123
GOODYEAR PD

Subject: TEST,TEST


DOB: 01/01/0101

Sex: M Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	06:13:47
Air Blank	0.000	06:14:15

*Improper Sample



A handwritten signature, possibly 'JR', is written in black ink. To the right of the signature, the number '#1123' is also handwritten in black ink.

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-000312

09/28/2006

06:14:42

Standard Lot#: 432401I

Last Changed By: J. VALLANDINGHAM #981

QAS: J. ROGERS #1123
GOODYEAR PD

Operator: J. ROGERS #1123
GOODYEAR PD

Subject: TEST,TEST

DOB: 01/01/0101

Sex: M Weight: 100

15 Minute Deprivation Period? Yes

Test	<u>g/210L</u>	Time
Air Blank	RFI*	06:15:17
Air Blank	0.000	06:15:46

*RFI Detect



Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-000312

09/28/2006

06:16:23

Standard Lot#: 432401I

Last Changed By: J. VALLANDINGHAM #981

QAS: J. ROGERS #1123
GOODYEAR PD

Operator: J. ROGERS #1123
GOODYEAR PD

Subject: TEST,TEST

DOB: 01/01/0101

Sex: M

Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	06:17:15
Diagnostic Test	Pass	06:17:42
Air Blank	0.000	06:18:09
0.100 Cal Check	0.098	06:18:28
Air Blank	0.000	06:18:58
Subject Test	0.017*	06:22:06
Air Blank	0.000	06:22:36
Air Blank	0.000	06:23:03
Subject Test	XXX**	06:23:38
Air Blank	0.000	06:24:11

*Deficient Sample

**Invalid Sample

Begin new deprivation period



#1123

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. ROGERS AGENCY GOODYEAR P.D.
DATE 9-28-06 TIME 0625
INTOXILYZER SERIAL # 80-000312 LOCATION GOODYEAR P.D.

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0.098 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

LOT 4324011
EXP 11/23/06

SIGNATURE [Signature] #1123

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number: 80-000312

09/28/2006

06:25:09


Standard Lot#: 432401I

Last Changed By: J. VALLANDINGHAM #981

QAS: J. ROGERS #1123

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	06:25:36
Cal Check	0.098	06:25:55
Air Blank	0.000	06:26:25



Operator's Signature