

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUELMANT 1144 AGENCY GOODYEAR PD  
DATE 9/3/2010 TIME 1812  
INTOXILYZER SERIAL # 80-000712 LOCATION GPA VENEZIA

- (R) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (M) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (A) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (S) 4. Air blank completed.
- (S) 5. Calibration check completed. Test results 0. 101 AC.
- (S) 6. Air blank completed.
- (S) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE  1144

INTOXILYZER 8000

Location: GOODYEAR PD VENIDA

Serial Number: 80-000312

Core Version: 8105.48

09/03/2010

18:12:28

Standard Lot#: 672733

Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144

GOODYEAR PD

Operator: M. GUILBAULT #1144

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	18:13:03
0.100 Cal Check	0.101	18:13:24
Air Blank	0.000	18:13:54

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. GUSTAFSON 1144 AGENCY GOOD YEAR PD  
DATE 9/3/10 TIME 1814  
INTOXILYZER SERIAL # 80-000312 LOCATION GPA VENEZA

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- 1. Clock time check.
- 2. Date check.

OPERATIONAL TESTS

- 1. Alcohol-free subject test result 0. 000 AC.
- 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE

M. Gustafson 1144

INTOXILYZER 8000  
Location: GOODYEAR PD VENIDA  
Serial Number: 80-000312  
Core Version: 8105.48  
09/03/2010

18:14:18

Standard Lot#: 672733  
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Subject: TEST,TEST  
DOB: 01/01/1990  
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:15:15
Diagnostic Test	Pass	18:15:50
Air Blank	0.000	18:16:18
0.100 Cal Check	0.104	18:16:39
Air Blank	0.000	18:17:09
Subject Test	0.000	18:17:34
Air Blank	RFI*	18:17:41
Air Blank	0.000	18:18:09

\*RFI Detect

Not a Successfully Completed Test Sequence



INTOXILYZER 8000  
Location: GOODYEAR PD VENIDA  
Serial Number: 80-000312  
Core Version: 8195.48  
09/03/2010

18:18:12

Standard Lot#: 672733  
Last Changed By: M. GUILBAULT #1144

OAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Subject: TEST, TEST  
DOB: 01/01/1920  
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	18:18:55
Air Blank	0.000	18:19:23

\*Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000  
Location: GOODYEAR PD VENIDA  
Serial Number: 80-000312  
Core Version: 8105.48  
09/03/2010

18:19:34

Standard Lot#: 672733  
Last Changed By: M. GUILBAULT #1144  
QAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Subject: TEST, TEST  
DOB: 01/01/1990  
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:20:35
Diagnostic Test	Pass	18:21:10
Air Blank	0.000	18:21:37
0.100 Cal Check	0.100	18:21:58
Air Blank	0.000	18:22:28
Subject Test	0.000	18:22:41
Air Blank	ART*	18:22:49
Air Blank	0.000	18:23:18

\*Sequence Aborted

Not a Successfully Completed Test Sequence

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. LAUREL MANT. 1144 AGENCY GOODY EAR PD  
DATE 9/13/2010 TIME 1832  
INTOXILYZER SERIAL # 80-000312 LOCATION GPS VENDA

- (1) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 099 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

INTOXILYZER 8000

Location: GOODYEAR PD VENTURA

Serial Number: 80-000312

Core Version: 8105.48

09/03/2000

18:25:00

Standard Lot#: 670723

Last Changed By: M. GUILBAULT #1141

QAS: M. GUILBAULT #1141  
GOODYEAR PD

Operator: M. GUILBAULT #1141  
GOODYEAR PD

Subject: TEST, TEST

DOB: 01/01/1990

Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	18:24:37
Diagnostic Test	Pass	18:25:12
Air Blank	0.000	18:28:19
0.100 Cal Check	0.101	18:36:00
Air Blank	0.000	18:36:30
Subject Test	0.072*	18:38:57
Air Blank	0.000	18:39:03
Air Blank	0.000	18:39:56
Subject Test	XXX**	18:31:03
Air Blank	0.000	18:31:37

\*Insufficient sample  
\*\*Invalid sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: GOODYEAR PD VERIDA

Serial Number: 30-000112

Firmware Version: #105.40

09/03/2010

12:32:01

Standard Lot#: 672733

Last Changed By: M. GUILBAULT #1144

LAB: M. GUILBAULT #1144

GOODYEAR PD

Operator: M. GUILBAULT #1144

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	12:32:31
0.100 Cal Check	0.999	12:32:52
Air Blank	0.000	12:33:21

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST W.C. Will #974 AGENCY Goodyear Police Department  
DATE 9/9/2010 TIME 0922  
INTOXILYZER SERIAL # 80-000312 LOCATION Goodyear Police Department

- W 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution AC.  
OR  
Pour a standard alcohol concentration solution 100 AC into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- W 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- W 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- W 4. Air blank completed.
- W 5. Calibration check completed. Test results 0 100 AC.
- W 6. Air blank completed.
- W 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J.P. Will

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GURBAJ 1144 AGENCY GOODHEAR PD  
DATE 9/18/10 TIME 1851  
INTOXILYZER SERIAL # 80-000312 LOCATION 190 N LITCHFIELD RD

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution -100 AC.  
OR

Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$

2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"  
 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.  
 4. Air blank completed.  
 5. Calibration check completed. Test results 0. 102 AC.  
 6. Air blank completed.  
 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

INTOXILYZER 8000  
Location: GPD 190 N LITCHFIELD  
Serial Number: 80-000312  
Core Version: 8105.48  
09/18/2010

18:51:33

-----  
Standard Lot#: 713873  
Last Changed By: M. GUILBAULT #1144  
-----

QAS: M. GUILBAULT #1144  
GOODYEAR PD  
-----

Operator: M. GUILBAULT #1144  
GOODYEAR PD  
-----

Test	g/210L	Time
Air Blank	0.000	18:52:04
0.100 Cal Check	0.102	18:52:25
Air Blank	0.000	18:52:55

Preventative Maintenance  
Performed -- Timer Reset



EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GURLO, AWC AGENCY GOOD YEAR PD  
DATE 9/20/2010 1144 TIME 1321  
INTOXILYZER SERIAL # 80-000212 LOCATION VENIDA

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 101 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE \_\_\_\_\_

INTOXILYZER 8000

Location: 175 N 145TH AVE, GYE

Serial Number: 80-000312

Core Version: 8105.48

09/20/2010

13:21:34

-----  
Standard Lot#: 767153

Last Changed By: M. GUILBAULT #1144

-----  
QAS: M. GUILBAULT #1144

GOODYEAR PD

-----  
Operator: M. GUILBAULT #1144

GOODYEAR PD

-----  
Test                      g/210L                      Time  
-----  
Air Blank                      0.000                      13:25:30  
0.100 Cal Check                0.101                      13:25:25  
Air Blank                      0.000                      13:25:20

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GURBAULT AGENCY Good Year PD  
DATE 9/20/10 TIME 1236  
INTOXILYZER SERIAL # 80-000312 LOCATION 175 N 145TH AVE

- (\*) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (\*) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (\*) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (\*) 4. Air blank completed.
- (\*) 5. Calibration check completed. Test results 0. 101 AC.
- (\*) 6. Air blank completed.
- (\*) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE \_\_\_\_\_

INTOXILYZER 8000  
Location: 175 N 145TH AVE, GYO  
Serial Number: 80-000312  
Core Version: 8105.48  
09/20/2010

12:36:23

Standard Lot#: 713873  
Last Changed By: M. GUILBAULT #1144

DAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	12:36:53
0.100 Cal Check	0.101	12:37:14
Air Blank	0.000	12:37:43

\*

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Will AGENCY Goodwin PD  
DATE 9/18/2010 TIME 1803  
INTOXILYZER SERIAL # 80-000312 LOCATION Goodwin

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"  
(✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.  
(✓) 4. Air blank completed.  
(✓) 5. Calibration check completed. Test results 0. 101 AC.  
(✓) 6. Air blank completed.  
(✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000

Location: GOODYEAR PD VENIDA

Serial Number: 80-000312

Core Version: 8105.48

09/18/2010

18:03:19

-----  
Standard Lot#: 672733

Last Changed By: M. GUILBAULT #1144  
-----

QAS: W. WILL #974

GOODYEAR PD  
-----

Operator: W. WILL #974

GOODYEAR PD  
-----

Test	g/210L	Time
Air Blank	0.000	18:03:50
0.100 Cal Check	0.101	18:04:11
Air Blank	0.000	18:04:40

Preventative Maintenance  
Performed -- Timer Reset