

INTOXILYZER 8000  
Instrument Initialization  
10:22:00 04/18/2012

INTOXILYZER 8000  
Location: 16161 W YUMA ROAD  
Serial Number: 80-001858  
Core Version: 8105.48  
04/18/2012 10:45:15

-----  
Standard Lot#: 34011100A1  
Last Changed By: M. GUILBAULT #1144

-----  
OAS: M. GUILBAULT #1144  
GOODYEAR PD

-----  
Operator: M. GUILBAULT #1144  
GOODYEAR PD

-----  
Test                    g/210L                    Time  
-----  
Air Blank                0.000                    10:45:45  
0.100 Cal Check        0.101                    10:46:06  
Air Blank                0.000                    10:46:35

Preventative Maintenance  
Performed -- Timer Reset


EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUILBAULT 1144 AGENCY GOODYEAR PD  
DATE 4/18/2012 TIME 1045  
INTOXILYZER SERIAL # 80-001858 LOCATION 16161 W Yuma RA

- (~~A~~) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (~~A~~) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"  
(~~A~~) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.  
(A) 4. Air blank completed.  
(~~A~~) 5. Calibration check completed. Test results 0. 101 AC.  
(~~A~~) 6. Air blank completed.  
(~~A~~) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE  1144

INTOXILYZER 8000  
Location: 16161 W YUMA ROAD  
Serial Number: 80-001858  
Core Version: 8105.48  
04/18/2012

*Handwritten mark*  
10:59:36

Standard Lot#: 34011100A1  
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Subject: TEST,TEST  
DOB: 01/01/1990  
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	11:00:31
Diagnostic Test	Pass	11:01:08
Air Blank	0.000	11:01:35
0.100 Cal Check	0.100	11:01:56
Air Blank	0.000	11:02:26
Subject Test	0.487*	11:05:33
Air Blank	0.000	11:06:06
Air Blank	0.000	11:06:34
Subject Test	XXX**	11:07:01
Air Blank	0.000	11:07:35

\*Deficient Sample  
\*\*Invalid Sample  
Begin new deprivation period

INTOXILYZER 8000  
Location: 16161 W YUMA ROAD  
Serial Number: 80-001858  
Core Version: 8105.48  
04/18/2012

*Handwritten mark*  
10:58:06

Standard Lot#: 34011100A1  
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Subject: TEST,TEST  
DOB: 01/01/1990  
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	ABT*	10:58:40
Air Blank	0.000	10:59:08

\*Sequence Aborted

Not a Successfully  
Completed  
Test Sequence

Not a Successfully  
Completed  
Test Sequence

INTOXILYZER 8000  
Instrument Initialization  
10:22:00 04/18/2012

INTOXILYZER 8000  
Location: 16161 W YUMA ROAD  
Serial Number: 80-001858  
Core Version: 8105.48  
04/18/2012 10:45:15

Standard Lot#: 34011100A1  
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	10:45:45
0.100 Cal Check	0.101	10:46:06
Air Blank	0.000	10:46:35

Preventative Maintenance  
Performed -- Timer Reset

INTOXILYZER 8000  
Location: 16161 W YUMA ROAD  
Serial Number: 80-001858  
Core Version: 8105.48  
04/18/2012

11:08:06

Standard Lot#: 34011100A1  
Last Changed By: M. GUILBAULT #1144

OAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	11:08:35
0.100 Cal Check	0.099	11:08:56
Air Blank	0.000	11:09:26

Preventative Maintenance  
Performed -- Timer Reset

INTOXILYZER 8000  
Location: 16161 W YUMA ROAD  
Serial Number: 80-001858  
Core Version: 8105.48  
04/18/2012

10:53:04

Standard Lot#: 34011100A1  
Last Changed By: M. GUILBAULT #1144

OAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Subject: TEST, TEST  
DOB: 01/01/1990  
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:55:14
Diagnostic Test	Pass	10:55:51
Air Blank	0.000	10:56:19
0.100 Cal Check	0.100	10:56:40
Air Blank	0.000	10:57:09
Subject Test	0.000	10:57:38
Air Blank	IPS*	10:57:45
Air Blank	RFI**	10:57:52

\*Improper Sample  
\*\*RFI Detect

Not a Successfully  
Completed  
Test Sequence

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. GUDLAULT AGENCY GOODYEAR PD  
DATE 2/18/12 TIME 1053  
INTOXILYZER SERIAL # 80-001858 LOCATION 16161 W Yuma Rd

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.  
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.  
 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

SIGNATURE [Signature] 1144

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. Gudubault AGENCY GOODYEAR PD  
DATE 4/18/12 TIME 1053 1108\*  
INTOXILYZER SERIAL # 80-001858 LOCATION 16161 W YUMA RD

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .120 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"  
(✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.  
(✓) 4. Air blank completed.  
(✓) 5. Calibration check completed. Test results 0. 099 AC.  
(✓) 6. Air blank completed.  
(✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE \_\_\_\_\_