

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. G. GIBSON AGENCY GOODYEAR PD  
DATE 2/4/2010 TIME 0139  
INTOXILYZER SERIAL # 82-001858 LOCATION SPD

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR

Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$

- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0. 101 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number: 80-001858

02/05/2010

01:38:43

Standard Lot#: 562576

Last Changed By: M. GUILBAULT #1144

OAS: M. GUILBAULT #1144

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	01:39:10
Cal Check	0.101	01:39:29
Air Blank	0.000	01:39:58

  
Operator's Signature

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUILBAULT 1144 AGENCY GOODYEAR PD  
DATE 2/24/2010 TIME 1619  
INTOXILYZER SERIAL # 80-001858 LOCATION GOODYEAR PD

- (\*) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (\*) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"  
(\*) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.  
(\*) 4. Air blank completed.  
(\*) 5. Calibration check completed. Test results 0. 101 AC.  
(\*) 6. Air blank completed.  
(\*) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

INTOXILYZER 8000  
Location: GOODYEAR POLICE DEPT  
Serial Number: 80-001858  
02/24/2010 16:19:02

-----  
Standard Lot#: 562576  
Last Changed By: M. GUILBAULT #1144  
-----

DAS: M. GUILBAULT #1144  
GOODYEAR PD  
-----

Test	g/210L	Time
Air Blank	0.000	16:19:30
Cal Check	0.101	16:19:49
Air Blank	0.000	16:20:18


  
-----  
Operator's Signature

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. GUILBAULT AGENCY Goodyear PD  
DATE 2/24/2010 TIME 1622  
INTOXILYZER SERIAL # 80-001858 LOCATION Goodyear PD

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- 1. Clock time check.
- 2. Date check.


OPERATIONAL TESTS

- 1. Alcohol-free subject test result 0. 000 AC.
- 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- 4. Standard Calibration Check standard 0. .100 AC. Result 0. 101 AC.

Instrument is operating properly and accurately. YES X NO     

COMMENTS

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE  1144

INTOXILYZER 8000  
Location: GOODYEAR POLICE DEPT  
Serial Number: 80-001858  
02/24/2010

16:22:24

Standard Lot#: 562576  
Last Changed By: M. GUILBAULT #1144

OAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Subject: TEST, TEST  
DOB: 10/10/2000  
Sex: M Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	16:23:24
Diagnostic Test	Pass	16:23:53
Air Blank	0.000	16:24:20
0.100 Cal Check	0.101	16:24:39
Air Blank	0.000	16:25:09
Subject Test	0.000	16:25:26
Air Blank	0.000	16:25:53
Five Minute Wait	ABT*	
Air Blank	0.000	16:26:29

Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000  
Location: GOODYEAR POLICE DEPT  
Serial Number: 80-001858  
02/24/2010

16:27:00

Standard Lot#: 562576  
Last Changed By: M. GUILBAULT #1144

OAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Subject: TEST, TEST  
DOB: 10/10/2000  
Sex: M Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	16:27:31
Air Blank	RFI**	16:27:40

\*Improper Sample  
\*\*RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000  
Location: GOODYEAR POLICE DEPT  
Serial Number: 80-001858  
02/24/2010

16:28:15

Standard Lot#: 562576  
Last Changed By: M. GUILBAULT #1144

OAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Subject: TEST, TEST  
DOB: 10/10/2000  
Sex: M Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	16:29:03
Diagnostic Test	Pass	16:29:32
Air Blank	0.000	16:30:00
0.100 Cal Check	0.101	16:30:18
Air Blank	0.000	16:30:48
Subject Test	0.014*	16:33:56
Air Blank	0.000	16:34:28
Air Blank	0.000	16:34:56
Subject Test	XXX**	16:35:19
Air Blank	0.000	16:35:53

\*Deficient Sample  
\*\*Invalid Sample  
Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUEBART 1144 AGENCY GOODYEAR PD  
DATE 2/24/2010 TIME 1639  
INTOXILYZER SERIAL # 80-001858 LOCATION GOODYEAR PD

- (~~1~~) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (~~2~~) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"  
(~~3~~) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.  
(~~4~~) 4. Air blank completed.  
(~~5~~) 5. Calibration check completed. Test results 0. 101 AC.  
(~~6~~) 6. Air blank completed.  
(~~7~~) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

INTOXILYZER 8000  
Location: GOODYEAR POLICE DEPT  
Serial Number: 80-001858  
02/24/2010 16:39:42

Standard Lot#: 562576  
Last Changed By: M. GUILBAULT #1144

DAS: M. GUILBAULT #1144  
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	16:40:10
Cal Check	0.101	16:40:28
Air Blank	0.000	16:40:58

  
Operator's Signature