

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUILBAULT AGENCY GOODYEAR PD  
DATE 6/16/2010 TIME 1443  
INTOXILYZER SERIAL # 80-001858 LOCATION GOODYEAR PD

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 101 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE  1144

INTOXILYZER 8000  
Location: GOODYEAR PD  
Serial Number: 80-001858  
Core Version: 8105.48  
06/16/2010

14:43:21

-----  
Standard Lot#: 526576  
Last Changed By: M. GUILBAULT #1144

-----  
QAS: M. GUILBAULT #1144  
GOODYEAR PD

-----  
Operator: M. GUILBAULT #1144  
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	14:43:53
0.100 Cal Check	0.101	14:44:14
Air Blank	0.000	14:44:44

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. GLEUBAUER AGENCY Goodyear PD  
DATE 6/16/10 TIME 1444  
INTOXILYZER SERIAL # 80-001858 LOCATION Goodyear PA

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.  
 2. Date check.

OPERATIONAL TESTS


1. Alcohol-free subject test result 0.          AC.  
 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 101 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE



1164

INTOXILYZER 8000  
Location: GOODYEAR PD  
Serial Number: 80-001858  
Core Version: 8105.48  
06/16/2010

14:44:58

---

Standard Lot#: 526576  
Last Changed By: M. GUILBAULT #1144

---

QAS: M. GUILBAULT #1144  
GOODYEAR PD

---

Operator: M. GUILBAULT #1144  
GOODYEAR PD

---

Subject: TEST,TEST  
DOB: 01/01/2000  
Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	14:46:05
Diagnostic Test	Pass	14:46:40
Air Blank	0.000	14:47:07
0.100 Cal Check	0.101	14:47:28
Air Blank	0.000	14:47:58
Subject Test	0.000	14:48:12
Air Blank	0.000	14:48:39
Wait	ABT*	
Air Blank	0.000	14:49:13

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: GOODYEAR PD  
Serial Number: 80-001858  
Core Version: 8105.48  
06/16/2010

14:49:56

Standard Lot#: 526576  
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Subject: TEST,TEST  
DOB: 01/01/2000  
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	14:50:30
Air Blank	RFI**	14:50:37

\*Improper Sample  
\*\*RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000  
Location: GOODYEAR PD  
Serial Number: 80-001858  
Core Version: 8105.48  
06/16/2010

14:59:17

---

Standard Lot#: 526576  
Last Changed By: M. GUILBAULT #1144

---

QAS: M. GUILBAULT #1144  
GOODYEAR PD

---

Operator: M. GUILBAULT #1144  
GOODYEAR PD

---

Subject: TEST,TEST

DOB: 01/01/2000

Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	15:00:06
Diagnostic Test	Pass	15:00:41
Air Blank	0.000	15:01:09
0.100 Cal Check	0.101	15:01:29
Air Blank	0.000	15:01:59
Subject Test	XXX*	15:02:24
Air Blank	0.000	15:02:58

\*Invalid Sample  
Begin new deprivation period

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: GOODYEAR PD  
Serial Number: 80-001858  
Core Version: 8105.48  
06/16/2010

14:50:40

---

Standard Lot#: 526576  
Last Changed By: M. GUILBAULT #1144

---

QAS: M. GUILBAULT #1144  
GOODYEAR PD

---

Operator: M. GUILBAULT #1144  
GOODYEAR PD

---

Subject: TEST,TEST  
DOB: 01/01/2000  
Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	14:51:30
Diagnostic Test	Pass	14:52:05
Air Blank	0.000	14:52:32
0.100 Cal Check	0.101	14:52:53
Air Blank	0.000	14:53:22
Subject Test	0.142*	14:56:30
Air Blank	0.000	14:57:02
Air Blank	0.000	14:57:29
Subject Test	0.163	14:58:32
Air Blank	ABT**	14:58:45
Air Blank	0.000	14:59:14

\*Deficient Sample  
\*\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M-GUELBAULT AGENCY GOODYEAR PD  
DATE 6/16/10 TIME ~~1444~~ 1503  
INTOXILYZER SERIAL # 80-001858 LOCATION GOODYEAR PD

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 101 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

1444



INTOXILYZER 8000  
Location: GOODYEAR PD  
Serial Number: 80-001858  
Core Version: 8105.48  
06/16/2010

15:03:13

-----  
Standard Lot#: 526576  
Last Changed By: M. GUILBAULT #1144  
-----

QAS: M. GUILBAULT #1144  
GOODYEAR PD  
-----

Operator: M. GUILBAULT #1144  
GOODYEAR PD  
-----

Test	g/210L	Time
Air Blank	0.000	15:03:43
0.100 Cal Check	0.101	15:04:04
Air Blank	0.000	15:04:33

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST W. GUTMANT AGENCY GOODYEAR PD  
DATE 6/24/2010 TIME 0104  
INTOXILYZER SERIAL # 80-0 0858 LOCATION 6PA #2

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0. 101 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000  
Location: GOODYEAR PD  
Serial Number: 80-001858  
Core Version: 8105.48  
06/25/2010


01:04:33

-----  
Standard Lot#: 526576  
Last Changed By: M. GUILBAULT #1144  
-----

QAS: M. GUILBAULT #1144  
GOODYEAR PD  
-----

Operator: M. GUILBAULT #1144  
GOODYEAR PD  
-----

Test	g/210L	Time
Air Blank	0.000	01:05:19
0.100 Cal Check	0.101	01:05:40
Air Blank	0.000	01:06:10

  
Preventative Maintenance  
Performed -- Timer Reset