

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUESSAULT 1144 AGENCY GOODYEAR PD  
DATE 3/15/2011 TIME 0627  
INTOXILYZER SERIAL # 80-001858 LOCATION 175 N 145TH AVE

- (A) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- ~~2~~ 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- ~~3~~ 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- ~~4~~ 4. Air blank completed.
- ~~5~~ 5. Calibration check completed. Test results 0. 101 AC.
- ~~6~~ 6. Air blank completed.
- ~~7~~ 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

INTOXILYZER 8000  
Location: 175 N 145TH AVENUE  
Serial Number: 80-001858  
Core Version: 8105.48  
03/15/2011

06:27:12

-----  
Standard Lot#: 767153  
Last Changed By: M. GUILBAULT #1144

-----  
QAS: M. GUILBAULT #1144  
GOODYEAR PD

-----  
Operator: M. GUILBAULT #1144  
GOODYEAR PD  
-----

Test	g/210L	Time
Air Blank	0.000	06:27:54
0.100 Cal Check	0.101	06:28:15
Air Blank	0.000	06:28:44

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. GUEMBOLT 1144 AGENCY GOODYEAR PD  
DATE 3/15/2011 TIME 0630  
INTOXILYZER SERIAL # 80-001858 LOCATION 175 N 145<sup>TH</sup> AVE

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS


1. Clock time check.  
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.  
 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 101 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

SIGNATURE 

INTOXILYZER 8000  
Location: 175 N 145TH AVENUE  
Serial Number: 80-001858  
Core Version: 8105.48  
03/15/2011

06:30:26 *HA*

---

Standard Lot#: 767153  
Last Changed By: M. GUILBAULT #1144

---

QAS: M. GUILBAULT #1144  
GOODYEAR PD

---

Operator: M. GUILBAULT #1144  
GOODYEAR PD

---

Subject: TEST,TEST  
DOB: 01/01/1990  
Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	06:31:21
Diagnostic Test	Pass	06:31:56
Air Blank	0.000	06:32:23
0.100 Cal Check	0.101	06:32:44
Air Blank	0.000	06:33:14
Subject Test	0.000	06:33:34
Air Blank	0.000	06:34:01
Wait	ABT*	
Air Blank	0.000	06:34:43

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: 175 N 145TH AVENUE  
Serial Number: 80-001858  
Core Version: 8105.48  
03/15/2011

06:34:49

---

Standard Lot#: 767153  
Last Changed By: M. GUILBAULT #1144

---

QAS: M. GUILBAULT #1144  
GOODYEAR PD

---

Operator: M. GUILBAULT #1144  
GOODYEAR PD

---

Subject: TEST,TEST  
DOB: 01/01/1990  
Sex: M Weight: 200

---

15 Minute Deprivation Period? Yes

---

Test	g/210E	Time
Air Blank	IPS*	06:35:23
Air Blank	RFI**	06:35:30

\*Improper Sample  
\*\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: 175 N 145TH AVENUE  
Serial Number: 80-001858  
Core Version: 8105.48  
03/15/2011

06:35:40

Standard Lot#: 767153  
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144  
GOODYEAR PD

Operator: M. GUILBAULT #1144  
GOODYEAR PD

Subject: TEST,TEST  
DOB: 01/01/1990  
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	06:36:52
Diagnostic Test	Pass	06:37:27
Air Blank	0.000	06:37:54
0.100 Cal Check	0.101	06:38:15
Air Blank	0.000	06:38:44
Subject Test	0.142*	06:41:52
Air Blank	0.000	06:42:24
Air Blank	0.000	06:42:52
Subject Test	XXX**	06:43:20
Air Blank	0.000	06:43:53

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUERMAUT 1144 AGENCY Goodyear PD  
DATE 3/15/2011 TIME 0644  
INTOXILYZER SERIAL # 80-001858 LOCATION 175 N 145TH AVE

- (D) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (D) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (D) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (D) 4. Air blank completed.
- (D) 5. Calibration check completed. Test results 0. 101 AC.
- (D) 6. Air blank completed.
- (D) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000  
Location: 175 N 145TH AVENUE  
Serial Number: 80-001858  
Core Version: 8105.48  
03/15/2011

06:44:15

-----  
Standard Lot#: 767153  
Last Changed By: M. GUILBAULT #1144

-----  
QAS: M. GUILBAULT #1144  
GOODYEAR PD

-----  
Operator: M. GUILBAULT #1144  
GOODYEAR PD  
-----

Test	g/210L	Time
Air Blank	0.000	06:44:44
0.100 Cal Check	0.101	06:45:05
Air Blank	0.000	06:45:34

Preventative Maintenance  
Performed -- Timer Reset