

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. COVELLBAULT 1144 AGENCY Goodyear PD
DATE 5/7/2010 TIME 1337
INTOXILYZER SERIAL # 80-001858 LOCATION 6PD

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0, 101 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
05/07/2010

17:37:46

Standard Lot#: 526576
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	17:38:15
0.100 Cal Check	0.101	17:38:36
Air Blank	0.000	17:39:06

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GIBLAVANT AGENCY GOODYEAR PD
DATE 5/7/2010 TIME 1650
INTOXILYZER SERIAL # 80-001858 LOCATION 6PD

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0.101 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

INVENTORY 8000
Location: GOODYEAR PD
Serial Number: E9-001328
Date Version: 0105.45
07/07/2018

16:50:00

Standard Code: 525036
Last Changed By: N. GOILBAULT #1144

QAS: N. GOILBAULT #1144
GOODYEAR PD

Operator: N. GOILBAULT #1144
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	16:51:10
0.100 Cal Check	0.101	16:51:31
Air Blank	0.000	16:51:54

Preventive Maintenance
Performed -- Timer Reset



INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
05/07/2010

16:42:56

Standard Lot#: 526576
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST, TEST
DOB: 01/01/1990
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	16:43:43
Diagnostic Test	Pass	16:44:18
Air Blank	0.000	16:44:46
0.100 Cal Check	0.101	16:45:07
Air Blank	0.000	16:45:36
Subject Test	0.071*	16:48:44
Air Blank	0.000	16:49:16
Air Blank	0.000	16:49:44
Subject Test	XXX**	16:49:58
Air Blank	0.000	16:50:31

~~*Deficient Sample~~

~~**Invalid Sample~~

Begin new deprivation period

A

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
05/07/2010

16:41:49

Standard Lot#: 526576
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST
DOB: 01/01/1990
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	16:42:26
Air Blank	RFI**	16:42:41

*Improper Sample ✓

**RFI Detect ✓

Not a Successfully Completed Test Sequence

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. GUELBAULT 1144 AGENCY GOODYEAR PD
DATE 5/7/2010 TIME 1637
INTOXILYZER SERIAL # 80-001853 LOCATION GPA

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 101 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE [Signature] 1144

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
05/07/2010

16:37:48

Standard Lot#: 526576
Last Changed By: M. GUILBAULT #1144

OAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST, TEST
DOB: 01/01/1990
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	16:38:36
Diagnostic Test	Pass	16:39:11
Air Blank	0.000	16:39:39
0.100 Cal Check	0.102	16:39:59
Air Blank	0.000	16:40:29
Subject Test	0.000	16:40:44
Air Blank	0.000	16:41:11
Wait	ABT*	
Air Blank	0.000	16:41:46

*Sequence Aborted

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUILBAULT 1144 AGENCY GOODYEAR PD
DATE 5/7/2010 TIME 1633
INTOXILYZER SERIAL # 80-001858 LOCATION GPA

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 102 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Cone Version: S105.48
95/07/2010

16:33:19

Standard Lot#: 526576
Last Changed By: M. GUILBAULT #1144

OAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	16:33:39
0.100 Cal Check	0.102	16:34:10
Air Blank	0.000	16:34:39

Preventative Maintenance
Performed -- Timer Reset



EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GLEBAUER 1144 AGENCY GOODYEAR PD
DATE 5/26/2010 TIME 0059
INTOXILYZER SERIAL # 80-001858 LOCATION GOODYEAR PD #2

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 100 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

1144

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
05/26/2010

00:59:18

Standard Lot#: 526576
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	00:59:51
0.100 Cal Check	0.102	01:00:12
Air Blank	0.000	01:00:41

Preventative Maintenance
Performed -- Timer Reset





**CERTIFICATE OF ANALYSIS
EBS - ETHANOL BREATH STANDARD**

CMI
316 EAST NINTH STREET
OWENSBORO, KY 42303

INVOICE#: 37420686
PO#: 374475
CMI ITEM #: 340135
DATE: 26 May 2010

ANALYTICAL ACCURACY: +/-0.002 BrAC or +/-2% whichever is greater.
PRODUCT EXPIRATION: 21 May 2012

CALGAZ LOT#: 767153
ETHANOL IN NITROGEN

COMPONENT	PPM	(BrAC)
ETHANOL	260.5	(0.100)
NITROGEN	BAL	

REFERENCE STANDARD	CYLINDER	CONCENTRATION PPM
NIST TRACEABLE STANDARDS*	D247819	259.6

* CERTIFICATION TRACEABLE TO N.I.S.T. RGM ETHANOL STANDARDS

MANUFACTURED DATE: 21 May 2010

CALGAZ CYLINDER SIZE: 6D

APPROVED BY : 

"We certify that all the cylinders for the Lot numbers identified herein are manufactured and tested within the requirements of CFR 49 part 178.65 and that physical and chemical test reports are on file and copies will be furnished upon request."

CALGAZ, Div. of Air Liquide Advanced Technologies U.S. LLC
821 Chesapeake Drive, Cambridge, MD 21613-0149
Phone: (410)228-6400 Fax: (410)228-4251