

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GORBAULT 1144 AGENCY GOODYEAR PD
DATE 11/12/2010 TIME 1354
INTOXILYZER SERIAL # 80-601858 LOCATION GPA / VENIQA

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 101 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE _____

[Signature] 1144

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
11/12/2010

13:54:18

Standard Lot#: 767153
Last Changed By: M. GUILBAULT #1144

QAS: M. CUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	13:55:04
0.100 Cal Check	0.101	13:55:25
Air Blank	0.000	13:55:54



Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. G. GILBAULT 1144 AGENCY GOODYEAR POLICE DEPT.
DATE 11/17/2010 TIME 0941
INTOXILYZER SERIAL # 80-001858 LOCATION GPD/VENIAA

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution, 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 100 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE _____

M. G. Gilbault 1144

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
11/17/2010

09:41:58

Standard Lot#: 767153
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	09:43:10
0.100 Cal Check	0.100	09:43:31
Air Blank	0.000	09:44:00

GA
1144

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. GUILBAULT 1144 AGENCY GOODYEAR POLICE DEPT.
DATE 11/17/2010 TIME 0950
INTOXILYZER SERIAL # 80-001858 LOCATION GPD/VENEDA

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 100 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE  1144

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
11/17/2010

09:44:43

Standard Lot#: 767153
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST
DOB: 01/01/1990
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:45:39
Diagnostic Test	Pass	09:46:14
Air Blank	0.000	09:46:42
0.100 Cal Check	0.100	09:47:02
Air Blank	0.000	09:47:32
Subject Test	0.000	09:47:46
Air Blank	IPS*	09:47:53
Air Blank	0.000	09:48:21

*Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
11/17/2010

09:48:33

Standard Lot#: 767153
Last Changed By: M. GUILBAULT #1144
QAS: M. GUILBAULT #1144
GOODYEAR PD
Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST
DOB: 01/01/1990
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:49:27
Diagnostic Test	Pass	09:50:02
Air Blank	0.000	09:50:30
0.100 Cal Check	0.100	09:50:51
Air Blank	0.000	09:51:20
Subject Test	0.000	09:51:39
Air Blank	0.000	09:52:06
Wait	ABT*	
Air Blank	0.000	09:52:41

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
11/17/2010

09:52:56

Standard Lot#: 767153
Last Changed By: M. GUILBAULT #1144

OAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST

DOB: 01/01/1990

Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:54:22
Diagnostic Test	Pass	09:54:57
Air Blank	0.000	09:55:25
Reference	RFI*	09:55:29
Air Blank	0.000	09:55:58

*RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
11/17/2010

09:56:03

Standard Lot#: 767153
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST
DOB: 01/01/1990
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:56:59
Diagnostic Test	Pass	09:57:34
Air Blank	0.000	09:58:02
0.100 Cal Check	0.100	09:58:22
Air Blank	0.000	09:58:52
Subject Test	0.033*	10:01:59
Air Blank	0.000	10:02:31
Air Blank	0.000	10:02:58
Subject Test	XXX**	10:03:15
Air Blank	0.000	10:03:49

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUILBAULT 1144 AGENCY GOOD YEAR POLICE DEPT.
DATE 11/17/2010 TIME 1004
INTOXILYZER SERIAL # 80-001858 LOCATION GPO/VENZAA

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 100 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE _____

M. Guilbault 1144

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
11/17/2010

10:04:15

Standard Lot#: 707153
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	10:04:43
0.100 Cal Check	0.100	10:05:00
Air Blank	0.000	10:05:35

Preventative Maintenance
Performed -- Timer Reset