

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST JEFF ROGERS AGENCY GOODYEAR P.D.
DATE 10-26-06 TIME 0719

INTOXILYZER SERIAL # 80-001858 LOCATION GOODYEAR JUSTICE FACILITY

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR

Pour a standard alcohol concentration solution 100 AC into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$

- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0 100 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

LOT 536401I
EXP 1-3-08


SIGNATURE [Signature] #1123

INTOXILIZER 8011
Location: GOODYEAR JUSTICE FAC
Serial Number: 80-001858
10/26/2016 07:19:36

Standard Lot#: 5364011
Last Changed By: K. SLONEKER H6559

DRS: J. ROGERS #1123
GOODYEAR PD

Test	9/210L	Time
Air Blank	0.000	07:20:05
Cal Check	0.100	07:20:24
Air Blank	0.000	07:20:54


Operator's Signature