

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. ROGERS AGENCY GOODYEAR P.D.
DATE 9-1-06 / 9-2-06 TIME 1746 / 0246
INTOXILYZER SERIAL # 80-001858 LOCATION GOODYEAR JUSTICE FACILITY

- TR* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution _____ AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- TR* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- TR* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- TR* 4. Air blank completed.
- TR* 5. Calibration check completed. Test results 0.111 / .100 AC.
- TR* 6. Air blank completed.
- TR* 7. Remove printed record. Attach the record to the completed checklist.


SIGNATURE *J. Rogers* #1123

INTOXILYZER 8000
Location: GOODYEAR JUSTICE FAC
Serial Number: 80-001858
09/01/2016 17:46:35

Standard Lot#: 5364011
Last Changed By: M. SLOVNER #6559

OAS: J. ROGERS #1123
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	17:46:33
Cal Check	0.100	17:46:52
Air Blank	0.000	17:47:22

 #1123


Operator's Signature

INTOXILYZER 8000
Location: GOODYEAR JUSTICE FAC
Serial Number: 80-001858
09/02/2016 02:46:18

Standard Lot#: 5364011
Last Changed By: M. SLOVNER #6559

OAS: J. ROGERS #1123
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	02:46:45
Cal Check	0.100	02:47:04
Air Blank	0.000	02:47:33

 #1123

Operator's Signature

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. ROGERS AGENCY GOODYEAR P.D.
DATE 9-29-06 TIME 1542
INTOXILYZER SERIAL # 80-001858 LOCATION GOODYEAR JUSTICE FACILITY

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution. 100 AC.
OR

Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$

2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
 4. Air blank completed.
 5. Calibration check completed. Test results 0. 101 AC. LOT 5364011
EXP 113/08
 6. Air blank completed.
 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] #1123

INTOXILYZER 3300
Location: GOODYER JUSTICE FRC
Serial Number: 90-001858
09/29/2006 15:42:48

Standard Lot#: 5364011
Last Changed By: M. SLOANEK #6559

DPS: J. ROGERS #1123
GOODYER PD

Test	9/21/06	Time
Rif Blank	0.000	15:43:15
Cal Check	0.101	15:43:34
Rif Blank	0.000	15:44:03



Operator Signature #1123

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST J. ROGERS AGENCY GOODYEAR P.D.
DATE 9-29-06 TIME 1546
INTOXILYZER SERIAL # 80-001858 LOCATION GOODYEAR JUSTICE FACILITY

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

LOT 5364011
EXP 1/3/08

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 102 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE [Signature] #1123

INTOXILYZER 8000
Location: GOODYEAR JUSTICE FAC
Serial Number: 80-001858
09/29/2006 16:33:54

Standard Lot#: 5364011
Last Changed By: M. SLOVEKER #6559

Operator: J. ROGERS #1123
GOODYEAR PD

Subject: TEST, TEST
DOB: 01/01/0101
Sex: M Weight: 100

15 Minute Depuration Period? Yes

Test g/210L Time

Air Blank 0.000 16:34:41
Diagnostic Test Pass 16:35:13
Air Blank 0.019 16:35:37
0.100 Cal Check 0.100 16:35:56
Air Blank 1.000 16:36:26
Subject Test 0.123* 16:36:34
Air Blank 0.000 16:40:06
Air Blank 0.000 16:40:33
Subject Test 0.000 16:41:28
Air Blank 0.000 16:42:00

*Deficient Sample
**Invalid Sample
Begin new depuration period

Not a Successfully
Completed
Test Sequence

INTOXILYZER 8000
Location: GOODYEAR JUSTICE FAC
Serial Number: 80-001858
09/29/2006 15:46:12

Standard Lot#: 5364011
Last Changed By: M. SLOVEKER #6559

Operator: J. ROGERS #1123
GOODYEAR PD

Subject: TEST, TEST
DOB: 01/01/0101
Sex: M Weight: 100

15 Minute Depuration Period? Yes

Test g/210L Time

Air Blank 0.000 15:47:07
Diagnostic Test Pass 15:47:36
Air Blank 0.000 15:48:03
0.100 Cal Check 0.101 15:48:22
Air Blank 0.000 15:48:52
Subject Test 0.000 15:49:31
Air Blank 0.000 15:49:59

*RFI Detect

Not a Successfully
Completed
Test Sequence

INTOXILYZER 8000
Location: GOODYEAR JUSTICE FAC
Serial Number: 80-001858
09/29/2006 16:25:24

Standard Lot#: 5364011
Last Changed By: M. SLOVEKER #6559

Operator: J. ROGERS #1123
GOODYEAR PD

Subject: TEST, TEST
DOB: 01/01/0101
Sex: M Weight: 100

15 Minute Depuration Period? Yes

Test g/210L Time

Air Blank 0.000 16:26:12
Diagnostic Test Pass 16:26:41
Air Blank 0.000 16:27:09
0.100 Cal Check 0.100 16:27:27
Air Blank 0.000 16:27:57
Subject Test 0.000 16:28:15
Air Blank 0.000 16:28:42
Five Minute Test 0.000 16:29:18

*Sequence Aborted

Not a Successfully
Completed
Test Sequence

INTOXILYZER 8000
Location: GOODYEAR JUSTICE FAC
Serial Number: 80-001858
09/29/2006 16:32:27

Standard Lot#: 5364011
Last Changed By: M. SLOVEKER #6559

Operator: J. ROGERS #1123
GOODYEAR PD

Subject: TEST, TEST
DOB: 01/01/0101
Sex: M Weight: 100

15 Minute Depuration Period? Yes

Test g/210L Time

Air Blank 0.000 16:33:04
Diagnostic Test Pass 16:33:33
Air Blank 0.000 16:33:33

*Improper Sample

Not a Successfully
Completed
Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. ROGERS AGENCY GOODYEAR P.D.

DATE 9-29-06 TIME 1643

INTOXILYZER SERIAL # 80-001858 LOCATION GOODYEAR JUSTICE FACILITY

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR

Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$

2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 102 AC. LOT 5364011
6. Air blank completed. EXP 1/3/08
7. Remove printed record. Attach the record to the completed checklist.


SIGNATURE [Signature] #1123

INITIALIZER 8011
Location: GOODYEAR JUSTICE PAC
Serial Number: 89-001858
09/25/2006 15:43:34

Standard LDR: 5364011
Last Changed By: M. SLOMKEK #6559

OPS: J. ROGERS #1123
GOODYEAR PG

Test	g/210L	Time
Air Blank	0.000	16:44:01
Cal Check	0.102	16:44:20
Air Blank	0.000	16:44:49


Operator's Signature