

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST K Miyazawa 1139 AGENCY GOODYEAR

DATE 5/1/2020 TIME 22:52

INTOXILYZER SERIAL # 80-000312 LOCATION 11 N 145TH AVE

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- 1. Clock time check.
- 2. Date check.

OPERATIONAL TESTS

- 1. Alcohol-free subject test result 0. 050 AC.
- 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- 4. Standard Calibration Check standard 0. 100 AC. Result 0. 096 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SIGNATURE [Signature] 1139

INTOXILYZER 8000  
Location: 11 N 145TH AVE  
Serial Number: 80-000312  
Core Version: 8105.48  
05/01/2020 23:07:52

Standard Lot#: 10818100A1  
Last Changed By: W. RYAN #1100

DAS: K. MIYAZATO #1139  
GOODYEAR PD

Operator: K. MIYAZATO #1139  
GOODYEAR PD

Subject: SOAP, TEST  
DOB: 04/02/1973  
Sex: M Weight: 203

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	23:08:40
Air Blank	0.000	23:09:08

\*Improper Sample

Not a Successfully  
Completed  
Test Sequence

INTOXILYZER 8000  
Location: 11 N 145TH AVE  
Serial Number: 80-000312  
Core Version: 8105.48  
05/01/2020 23:11:26

Standard Lot#: 10818100A1  
Last Changed By: W. RYAN #1100

DAS: K. MIYAZATO #1139  
GOODYEAR PD

Operator: K. MIYAZATO #1139  
GOODYEAR PD

Subject: SOAP, TEST  
DOB: 04/02/1973  
Sex: M Weight: 203

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	23:12:43
Diagnostic Test	Pass	23:13:20
Air Blank	0.000	23:13:47
0.100 Cal Check	0.096	23:14:08
Air Blank	0.000	23:14:38
Subject Test	0.026*	23:17:45
Air Blank	0.000	23:18:16
Air Blank	0.000	23:18:46
Subject Test	XXX**	23:19:13
Air Blank	0.000	23:19:46

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

Not a Successfully  
Completed  
Test Sequence

INTOXILYZER 8000  
Location: 11 N 145TH AVE  
Serial Number: 80-000312  
Core Version: 8105.48  
05/01/2020 23:00:13

Standard Lot#: 10818100A1  
Last Changed By: W. RYAN #1100

DAS: K. MIYAZATO #1139  
GOODYEAR PD

Operator: K. MIYAZATO #1139  
GOODYEAR PD

Subject: SCAP, TEST  
DOB: 04/02/1973  
Sex: M Weight: 203

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	23:01:33
Diagnostic Test	Pass	23:02:11
Air Blank	0.000	23:02:38
0.100 Cal Check	0.097	23:02:59
Air Blank	0.000	23:03:29
Subject Test	0.000	23:03:51
Air Blank	ABT*	23:04:05
Air Blank	0.000	23:04:33

\*Sequence Aborted

Not a Successfully  
Completed  
Test Sequence

INTOXILYZER 8000  
Location: 11 N 145TH AVE  
Serial Number: 80-000312  
Core Version: 8105.48  
05/01/2020 23:09:29

Standard Lot#: 10818100A1  
Last Changed By: W. RYAN #1100

DAS: K. MIYAZATO #1139  
GOODYEAR PD

Operator: K. MIYAZATO #1139  
GOODYEAR PD

Subject: SCAP, TEST  
DOB: 04/02/1973  
Sex: M Weight: 203

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	23:10:11
Air Blank	0.000	23:10:39

\*RFI Detect

Not a Successfully  
Completed  
Test Sequence

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST K Miyazato 1139 AGENCY Goodyear  
DATE 5/1/2020 TIME 23:20  
INTOXILYZER SERIAL # 86-000312 LOCATION 11 N 145TH AVE

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0.096 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1139

INTOXILYZER 8000  
Location: 11 N 145TH AVE  
Serial Number: 60-000312  
Core Version: 8103.46  
15/01/2020 23:20:23

Standard Lot#: 10818100A1  
Last Changed By: W. RYAN #1100

DAS: K. MIYAZATO #1139  
GOODYEAR PD

Operator: K. MIYAZATO #1139  
GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	23:21:54
0.100 Cal Check	0.096	23:21:15
Air Blank	0.000	23:21:45

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST K MLYNARD 1139 AGENCY GOODYEAR  
DATE 5/30/20 TIME 2117  
INTOXILYZER SERIAL # 80-000312 LOCATION 11 N 145TH AVE

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 096 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1139

INTOXILYZER 8000  
Location: 11 N 145TH AVE  
Serial Number: 80-000312  
Core Version: 8105.48  
05/30/2020 21:17:52

-----  
Standard Lot#: 10818100A1  
Last Changed By: W. RYAN #1100  
-----

QAS: K. MIYAZATO #1139  
GOODYEAR PD  
-----

Operator: K. MIYAZATO #1139  
GOODYEAR PD  
-----

Test	g/210L	Time
Air Blank	0.000	21:18:23
0.100 Cal Check	0.096	21:18:44
Air Blank	0.000	21:19:15

Preventative Maintenance  
Performed -- Timer Reset