

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUELBAUT 1144 AGENCY GOOSEYAN PD
DATE 9/3/2010 TIME 1518
INTOXILYZER SERIAL # 80-001858 LOCATION GPD #2

- ~~A~~ 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- ~~A~~ 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- ~~A~~ 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- ~~A~~ 4. Air blank completed.
- ~~A~~ 5. Calibration check completed. Test results 0. 101 AC.
- ~~A~~ 6. Air blank completed.
- ~~A~~ 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE _____

M. Guelbaut 1144

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
09/03/2010

15:18:01 *A*

Standard Lot#: 526576
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

| Test | g/210L | Time |
|-----------------|--------|-------------------|
| Air Blank | 0.000 | 15:18:32 |
| 0.100 Cal Check | 0.101 | 15:18:53 <i>A</i> |
| Air Blank | 0.000 | 15:19:23 |

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. BURBANT 1144 AGENCY GOODYEAR PD

DATE 9/3/2010 TIME 1520

INTOXILYZER SERIAL # 80-001858 LOCATION GPD#2

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.
() 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.
() 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
() 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
() 4. Standard Calibration Check standard 0. 100 AC. Result 0. 095 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE  1144

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
09/03/2010

15:20:49

Standard Lot#: 526576
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST, TEST
DOB: 01/01/1990
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

| Test | g/210L | Time |
|-----------------|--------|----------|
| Air Blank | 0.000 | 15:21:45 |
| Diagnostic Test | Pass | 15:22:20 |
| Air Blank | 0.000 | 15:22:47 |
| 0.100 Cal Check | 0.101 | 15:23:08 |
| Air Blank | 0.000 | 15:23:38 |
| Subject Test | 0.000 | 15:23:51 |
| Air Blank | ABT* | 15:24:01 |
| Air Blank | 0.000 | 15:24:30 |

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
09/03/2010

15:24:35

Standard Lot#: 526576
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST
DOB: 01/01/1990
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

| Test | g/210L | Time |
|-----------|-----------------|----------|
| Air Blank | IPS* | 15:25:03 |
| Air Blank | RFI** | 15:25:12 |

*Improper Sample
**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
09/03/2010

15:25:24

Standard Lot#: 526576
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

Subject: TEST,TEST
DOB: 01/01/1990
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

| Test | g/210L | Time |
|-----------------|--------|----------|
| Air Blank | 0.000 | 15:26:11 |
| Diagnostic Test | Pass | 15:26:46 |
| Air Blank | 0.000 | 15:27:13 |
| 0.100 Cal Check | 0.101 | 15:27:34 |
| Air Blank | 0.000 | 15:28:04 |
| Subject Test | 0.222* | 15:31:11 |
| Air Blank | 0.000 | 15:31:43 |
| Air Blank | 0.000 | 15:32:11 |
| Subject Test | XXX** | 15:32:23 |
| Air Blank | 0.000 | 15:32:57 |

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST R. Guebara AGENCY Goodyear PD
DATE 9/3/2010 TIME 1534
INTOXILYZER SERIAL # 80-001858 LOCATION GPD #2

- (S) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution, 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (S) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (S) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (S) 4. Air blank completed.
- (S) 5. Calibration check completed. Test results 0. 095 AC.
- (S) 6. Air blank completed.
- (S) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

INTOXILYZER 8000

Location: GOODYEAR PD

Serial Number: 80-001858

Core Version: 8105.48

09/03/2010

15:34:53

Standard Lot#: 526576

Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144

GOODYEAR PD

Operator: M. GUILBAULT #1144

GOODYEAR PD

| Test | g/210L | Time |
|-----------------|--------|----------|
| Air Blank | 0.000 | 15:35:23 |
| 0.100 Cal Check | 0.095 | 15:35:44 |
| Air Blank | 0.000 | 15:36:13 |

Preventative Maintenance
Performed -- Timer Reset

INTOXILYZER 8000

Location: GOODYEAR PD

Serial Number: 80-001858

Core Version: 8105.48

09/09/2010

09:22:50

Standard Lot#: 526576

Last Changed By: M. GUILBAULT #1144

QAS: W. WILL #974

GOODYEAR PD

Operator: W. WILL #974

GOODYEAR PD

| Test | g/210L | Time |
|-----------------|--------|----------|
| Air Blank | 0.000 | 09:23:20 |
| 0.100 Cal Check | 0.102 | 09:23:41 |
| Air Blank | 0.000 | 09:24:11 |

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE


QA SPECIALIST M. GUTENBAIT 1144 AGENCY GOODYEAR PD
DATE 9/20/2010 TIME 1646
INTOXILYZER SERIAL # 80-001858 LOCATION GPA #2

- (*) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (*) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (*) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (*) 4. Air blank completed.
- (*) 5. Calibration check completed. Test results 0. 101 AC.
- (*) 6. Air blank completed.
- (*) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
09/20/2010

16:46:10



Standard Lot#: 526576
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

| Test | g/210L | Time |
|-----------------|--------|----------|
| Air Blank | 0.000 | 16:46:41 |
| 0.100 Cal Check | 0.101 | 16:47:02 |
| Air Blank | 0.000 | 16:47:31 |

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. Galtman 1144 AGENCY GOODYEAR PA
DATE 9/20/2010 TIME 1650
INTOXILYZER SERIAL # 80-001858 LOCATION GPS #2

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution, 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 101 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE M. Galtman 1144

INTOXILYZER 8000
Location: GOODYEAR PD
Serial Number: 80-001858
Core Version: 8105.48
09/20/2010

16:50:38

Standard Lot#: 767153
Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144
GOODYEAR PD

Operator: M. GUILBAULT #1144
GOODYEAR PD

| Test | g/210L | Time |
|-----------------|--------|----------|
| Air Blank | 0.000 | 16:51:07 |
| 0.100 Cal Check | 0.101 | 16:51:28 |
| Air Blank | 0.000 | 16:51:58 |

Preventative Maintenance
Performed -- Timer Reset