

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUILHAULT AGENCY GOODYEAR PD
DATE 9/9/09 TIME 0856
INTOXILYZER SERIAL # 80-001858 LOCATION GOODYEAR PD

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 101 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number: 80-001858

09/09/2009

08:55:47

Standard Lot#: 562576

Last Changed By: M. GUILBAULT #1144

OAS: M. GUILBAULT #1144

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	08:56:14
Cal Check	0.101	08:56:33
Air Blank	0.000	08:57:02



Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. GUELMAN AGENCY GOODYEAR PD
DATE 9/9/09 TIME 0922
INTOXILYZER SERIAL # 80-001853 LOCATION GOODYEAR PD

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE [Signature]

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-001858

09/09/2009

09:20:45

Standard Lot#: 562576

Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144

GOODYEAR PD

Operator: M. GUILBAULT #1144

GOODYEAR PD

Subject: SQAP, TEST

DOB: 01/01/2000

Sex: M

Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:22:10
Diagnostic Test	Pass	09:22:38
Air Blank	0.000	09:23:05
0.100 Cal Check	0.100	09:23:24
Air Blank	0.000	09:23:53
Subject Test	0.000	09:24:11 ✓
Air Blank	0.000	09:24:38 ✓
Five Minute Wait	ABT*	
Air Blank	0.000	09:25:26

*Sequence Aborted

Not a Successfully Completed Test Sequence ✓

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-001858

09/09/2009

09:27:26

Standard Lot#: 562576

Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144

GOODYEAR PD

Operator: M. GUILBAULT #1144

GOODYEAR PD

Subject: SQAP,TEST

DOB: 01/01/2000

Sex: M

Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:28:32
Diagnostic Test	Pass	09:29:00
Air Blank	IPS*	09:29:06
Air Blank	RFI**	09:29:30

*Improper Sample

**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number : 80-001858

09/09/2009

09:30:27

Standard Lot#: 562576

Last Changed By: M. GUILBAULT #1144

QAS: M. GUILBAULT #1144

GOODYEAR PD

Operator: M. GUILBAULT #1144

GOODYEAR PD

Subject: SQAP,TEST

DOB: 01/01/2000

Sex: M

Weight: 100

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	09:31:32
Diagnostic Test	Pass	09:31:59
Air Blank	0.000	09:32:26
0.100 Cal Check	0.100	09:32:45
Air Blank	0.000	09:33:15
Subject Test	0.358*	09:36:23 ✓
Air Blank	0.000	09:37:00
Air Blank	0.000	09:37:27 ✓
Subject Test	XXX**	09:38:07
Air Blank	0.000	09:38:46

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. GUEUBAULT AGENCY GOODYEAR PD
DATE 9/9/09 TIME 0940
INTOXILYZER SERIAL # 80-001898 LOCATION GOODYEAR PD

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 099 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature] 1144

INTOXILYZER 8000

Location: GOODYEAR POLICE DEPT

Serial Number: 80-001858

09/09/2009

09:40:09

Standard Lot#: 562576

Last Changed By: M. GULLBAULT #1144

QAS: M. GULLBAULT #1144

GOODYEAR PD

Test	g/210L	Time
Air Blank	0.000	09:40:36
Cal Check	0.099	09:40:55
Air Blank	0.000	09:41:25


Operator's Signature