

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST W. T. RYAN #1100 AGENCY GOODYEAR POLICE DEPT  
DATE 5-29-17 TIME 0133  
INTOXILYZER SERIAL # 80-001854 LOCATION 11 N 145TH AVE GOODYEAR AZ

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- ( ) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 099 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE W. T. Ryan #1100

GOODMAN FORM 100  
0133

W. T. Page #100  
2-28-10

11 N. 14TH AVE. GOODMAN NE  
0010

80-001820

000

W. T. Page #100

Handwritten vertical text on the right margin, possibly a date or reference number.

INTOXILYZER 8000  
Location: 11 S 145TH AVE  
Serial Number: 80-001858  
Core Version: 8105.48  
05/29/2017 01:33:02

-----  
Standard Lot#: 28515100A1  
Last Changed By: W. RYAN #1100  
-----

QAS: W. RYAN #1100  
GOODYEAR PD  
-----

Operator: W. RYAN #1100  
GOODYEAR PD  
-----

Test	g/210L	Time
Air Blank	0.000	01:33:38
0.100 Cal Check	0.099	01:33:59
Air Blank	0.000	01:34:29

Preventative Maintenance  
Performed -- Timer Reset

*W. Ryan*